

cytomin^e and BIAflows for data and computer scientists

May 11th 2021

Presenters : Raphaël Marée & Sébastien Tosi




Helpers : Renaud Hoyoux, Ulysse Rubens, Romain Mormont



Last week: Cytomine-webUI for collaborative annotation


2

Video recording: <https://www.youtube.com/watch?v=AFcILoaaekw>





Single-Page Application, Written in VueJS, loaded in browser:
<https://github.com/cytomine/Cytomine-Web-UI>

Webinar PART 1



Webinar PART 2



Cytomine

Overview of existing Apps

Architecture

Data models

API

API Clients

Java/Javascript/Python

Creating Reproducible Apps

Descriptor

Containers

BIAFLOWS

Benchmarking with demo

Running algorithms from Cytomine Web-UI:

1. Configure “trusted sources” to fetch “Apps” from a source

The screenshot shows the Cytomine Web-UI interface. On the left, a modal dialog titled "Update trusted source" is open, showing fields for "Source code provider username" (cytome), "Environment provider username" (cytome), "Prefix" (S_), "Source code provider" (GitHub), and "Environment provider" (docker). On the right, the main dashboard displays the "Cytomine" project page, which includes sections for "Pinned repositories" (Cytomine-bootstrap, Cytomine-core) and a list of other repositories like Cytomine-Web-UI and S_Segment-ML-MultiOutputET-Pred-BI.

Configuration to be done at the level of your Cytomine server (by admin e.g. core facility manager, PI)

Running algorithms from Cytomine Web-UI:


2. Enable “Apps” (“Algorithms”) in your project configuration

General	Members	Custom UI	Algorithms	Image filters
<input type="text"/> Search				
Name		Version	Runnable	Status
Stats-TermArea (v1.1.1)		Last release	✓ Yes	Enabled
CellDetect_Stardist_HE_ROI (v1.0.3)		Last release	✓ Yes	Enabled
Segment-CV-AdaptThres-Sample (v1.3.0)		Last release	✓ Yes	Enabled
Segment-CV-AdaptThres-Sample (v1.2.1)		Deprecated	✓ Yes	Enabled
Segment-CV-AdaptThres-Sample (v1.2.0)		Deprecated	✓ Yes	Enabled
BckgrdSeg-ML-RandSubwET-Pred-BI		No information	✗ No	Enabled
Segment-CV-AdaptThres-Object-BI (v1.0.3)		Last release	✓ Yes	Enabled
Segment-CV-AdaptThres-Sample (v1.1.0)		Deprecated	✓ Yes	Enabled
TissueSegment_Model_Predict		No information	✗ No	Enabled
TissueSegment_Model_Builder		No information	✗ No	Enabled
TissueDetect		No information	✗ No	Enabled
S_Zebrafish_Head_Operculum_UNet_Segmentation (v3.0.4)		Last release	✓ Yes	Disabled

Configuration to be done at the project level (by project manager e.g. PI or researcher)

Example: Tissue sample detection (no AI, basic image processing)⁶

- Code: https://github.com/cytomine/S_Segment-CV-AdaptThres-Sample
- Input: List of images, ontology term to predict
- Output: Annotations (contours) of tissue samples with ontology term



Launch new analysis


Algorithm Segment-CV-AdaptThres-Sample (v1.3.0)

Name	Value
Term to predict	Select options
Pre-filled parameters Hide	
Maximum image size	2048
Threshold block size	951
Threshold constant	5
Number of erosions	3
Number of dilations	3
Level of log	INFO

[Cancel](#) [Launch new analysis](#)

Example: Tissue sample detection (no AI, basic image processing)⁷


- Code: https://github.com/cytomine/S_Segment-CV-AdaptThres-Sample
- Input: List of images, ontology term to predict
- Output: Annotations (contours) of tissue samples with ontology term



Example: Object detection in hyperspectral (2D+c) images

8

- Code: https://github.com/Cytomine-ULiege/S_Segment-CV-Object-Projection
- Input: List of images
- Output: Annotations (contours) of “objects”



Launch new analysis

Algorithm Segment-CV-Object-Projection (v1.0)


Name	Value
Images to process	G1_7_top.czi
Optional parameters Hide	
Term to predict	ROI
Pre-filled parameters Hide	
Projection	max
Thresholding filter	otsu
Tile size	1024
Tile overlap	32
Minimum Object Area	100
Slices to use for annotations	median

Cancel Launch new analysis

Example: Object detection in hyperspectral (2D+c) images

9


- Code: https://github.com/Cytomine-ULiege/S_Segment-CV-Object-Projection
- Input: List of images
- Output: Annotations (contours) of “objects”



(Using Maximum Intensity Projection & Otsu)

Example: Cell segmentation (pre-trained StarDist H&E model)

- Code: https://github.com/cytomine/S_CellDetect_Stardist_HE_ROI
- Inputs: List of WSIs & ontology term for regions of interest, Stardist parameters
- Outputs: in each ROI, cell segmentations with ontology term




(based on Schmidt & Weigert, MICCAI 2018;
<https://github.com/stardist/stardist>)

Example: Tissue area delineation: training (Random forests)

11

- Code: https://github.com/cytomine/S_Segment-ML-MultiOutputET-Train
- Inputs: images & ontology terms of positive/negative annotations, RF parameters (nb of trees, subwindow sizes,...)
- Outputs: a random forest-based segmentation model (saved in Cytomine database as “AttachedFile”)




(Marée et al., ISBI 2014)

Example: Tissue area delineation: prediction (Random forests)

12


- Code: https://github.com/cytomine/S_Segment-ML-MultiOutputET-Pred-BI
- Inputs: list of images, a software identifier which generated a model, ontology term
- Outputs: in each image, segmented positive regions with ontology term



(Marée et al., ISBI 2014)

Example: Foreground segmentation (U-Net)


- Code: https://github.com/Cytomine-ULiege/S_Segment-ML-UNet-Binary-Train
https://github.com/Cytomine-ULiege/S_Segment-ML-UNet-Binary-Pred
- Inputs: ontology terms for foreground annotations, U-Net parameters (nb. epochs, learning rate, weight decay, ...)
- Outputs: a U-Net segmentation model (saved in Cytomine database as “[AttachedFile](#)”)



(Ronneberger, et al., ISBI 2015)



Example: Zebrafish Head/Operculum: prediction (pre-trained U-Net)¹⁴

- Code:
https://github.com/navdeepkaushish/S_Zebrafish_Head_Operculum_UNet_Segmentation
- Inputs: list of images, ontology terms for head & operculum
- Outputs: in each image, segmented head and operculum regions with corresponding terms



Example: Landmark detection : training (Random forests)




- Code: https://github.com/Cytomine-ULiege/S_LandmarkDetect-ML-MSET-Train
- Inputs: images & ontology term of manual annotations of anatomical landmarks
- Outputs: a random forest based detection model for each landmark (saved in Cytomine database as “AttachedFile”)



(Vandaele et al., Scientific Reports, 2019)

Example: Landmark detection : prediction (Random forests)

- Code: https://github.com/Cytomine-ULiege/S_LandmarkDetect-ML-MSET-Pred
- Inputs: images & ontology terms of manual annotations of anatomical landmarks
- Outputs: predicted annotation points for landmarks



Example: Statistics exporter

- Code: https://github.com/Cytomine-ULiege/S_Stats-TermArea
- Inputs: List of images, list of ontology terms, user/software ids
- Outputs: Statistics (counts, area, ratio, ...)

Stats-TermArea (v1.1.2) 1 rmaree Mar 26, 2021 4:17 PM



Status comment	Finished										
Execution duration	a few seconds										
Parameters	Show										
Execution log	Show										
Tags	No tag Add										
Files	<table border="1"> <thead> <tr> <th>Filename</th> <th>Comment</th> <th>Size</th> <th></th> <th></th> </tr> </thead> <tbody> <tr> <td>stats.csv</td> <td>output</td> <td>759 B</td> <td>View</td> <td>Download</td> </tr> </tbody> </table>	Filename	Comment	Size			stats.csv	output	759 B	View	Download
Filename	Comment	Size									
stats.csv	output	759 B	View	Download							



A	B	C	D
1 Area statistics			
2 Image	Nucleus	ROI	Total
3 CMU-1.tif	362989	362989	725978
4 Total	362989	362989	725978
5			
6			
7 Number statistics			
8 Image	Nucleus	ROI	Total
9 CMU-1.tif	4	4	8
10 Total	4	4	8
11			
12			
13 Data ratio			
14 Image	Nucleus	ROI	Total
15 CMU-1.tif	0.5	0.5	1
16			
17			
18 Details			
19 *****			
20 Image 1			
21 CMU-1.tif			
22 *****			
23 #####			
24 Nucleus			
25 Created	Area		
26 2021.03.19 10:44:23	27689		
27 2021.03.19 10:43:54	0		
28 2021.03.19 09:50:49	0		
29 2021.03.11 09:57:26	335300		
30			
31 Number of annotations	4		
32 Average area	90747		
33 Total area	362989		
34			
35			
36 #####			
37 ROI			
38 Created	Area		
39 2021.03.19 10:44:23	27689		
40 2021.03.19 10:43:54	0		
41 2021.03.19 09:50:49	0		
42 2021.03.11 09:57:26	335300		
43			
44 Number of annotations	4		
45 Average area	90747		
46 Total area	362989		

Results are not satisfactory ? → Reviewing of algorithms predictions¹⁸

- Results (e.g. cell detections, tissue delineation) can be reviewed/proofread (original & validated annotations are stored in Cytomine database)




Results are not satisfactory ?

19

→ Develop and integrate into Cytomine your own (better) algorithms !




Beyond the scene: Cytomine architecture (simplified)



The data stored and managed by a Cytomine server can be obtained through specific URLs

Beyond the scene: Cytomine architecture (simplified)




*All the endpoints are secured.
You need to be authenticated
so that the Cytomine server
accept to treat your request.*

Data sources:

- SQL & NoSQL Databases (Docker volumes)
- Image files (Filesystem)

Beyond the scene: Cytomine architecture & libraries



Leveraging **Docker** technology, you can **install** the whole Cytomine platform on a personal laptop (isolated so without collaborative features), on a virtual machine or a dedicated server in your lab or at large scale (e.g. core facility or MOOC for thousands of users), by hosting containers on different machines to increase performance and stability (horizontal scaling).


<https://doc.cytomine.org/admin-guide/>

“Toy” demo server:
<https://demo.cytomine.be>
username: jsnow
password: jsnow

To request a specific demo server access: info@cytomine.org


Beyond the scene: Cytomine data models (simplified)

- User: Me and You and Everyone We Know
- Permission:
 - Manage access rights
- Project:
 - Shared area between users
 - Group images & annotations
- Ontology:
 - Hierarchical list of terms to be associated to annotations
- Annotation:
 - Geometry + term within an image
- Image:
 - Reference to underlying image file (2D+c+z+t)
- Storage:
 - Virtual disk space
- Metadata:
 - AttachedFile: file attached as a metadata
 - Description: rich text description
 - Property: key-value pairs
 - Tags
- Software:
 - Runnable tools (e.g AI algorithm) from Cytomine



Beyond the scene: Cytomine annotation data model (simplified)

- 3 types
 - UserAnnotation: done by a human user
 - 1 annot layer per user in an image
 - AlgoAnnotation: done by a software execution (UserJob)
 - 1 annot layer per execution in an image
 - ReviewedAnnotation:
 - User or algo annotations
 - Reviewed by an “expert”
 - 0/1 review layer in an image
- Geometry format: WKT
https://en.wikipedia.org/wiki/Well-known_text_representation_of_geometry
- Stored in a spatial database (PostGIS)




- + AnnotationLinks to group annotations depicting same regions in multiple modalities: Under development (Cytomine ULiège R&D)

Cytomine Interoperability: HTTPS REST API

- API endpoints
 - Always require authentication
 - JSON

- Example: Project resource

- List projects I have access to
 - GET <https://demo.cytomine.be/api/project.json>
 - Support filters like minImages, minAnnots, ...
- Get a specific project (if I have rights)
 - GET <https://demo.cytomine.be/api/project/528050.json>
- Update project info (if I have rights)
 - PUT <https://demo.cytomine.be/api/project/80741.json> + JSON body
- Delete project (if I have rights)
 - DELETE <https://demo.cytomine.be/api/project/80741.json>
- Add new project (if I have rights)
 - POST <https://demo.cytomine.be/api/project.json> + JSON body



```
{
  class: "be.cytomine.project.Project",
  id: 80741,
  created: "1615273192548",
  updated: "1617106983480",
  deleted: null,
  name: "DEMO-VARIOUS",
  ontology: 80733,
  ontologyName: "DEMO-VARIOUS",
  discipline: null,
  blindMode: false,
  areImagesDownloadable: false,
  disciplineName: null,
  number_of_slides: 0,
  number_of_images: 8,
  number_of_annotations: 46,
  number_of_job_annotations: 4108,
  retrieval_projects: [ ],
  number_of_reviewed_annotations: 2,
  retrieval_disable: false,
  retrieval_all_ontology: true,
  is_closed: false,
  is_read_only: false,
  is_restricted: false,
  hide_users_layers: false,
  hide_admins_layers: false
}
```

- Same idea for other Cytomine resources: images, annotations, ontologies, users, ...
- **Anything that can be done from Cytomine-WebUI can be done using the HTTPS REST API ... and even more than that!**


Annotations in Cytomine REST API

- List annotations based on filters (image, user, ontology term,...):

<https://demo.cytomine.be/api/annotation/search.json?image=528132&user=701&term=528044>

- Get a specific annotation:

<https://demo.cytomine.be/api/annotation/143584607.json>




Annotation Geometry in WKT format




There are WKT parsers to convert WKT
from and to geometry objects in several
languages (e.g. Python: Shapely)

```
{
  "class": "be.cytomine.ontology.UserAnnotation",
  "id": 143584607,
  "created": "1575556669131",
  "updated": "1575556707121",
  "deleted": null,
  "slice": 146632788,
  "project": "528050",
  "image": 528132,
  "user": 701,
  "location": "POLYGON ((23569.544251797495 12583.846593241735, 23760.89049605907 12583.846593241735, 23827.23647439271 12572.308144907995, 23871.46743142396 12551.154313426061, 23916.659692123966.659615872202 12500.192742686315, 24008.96727883607 12473.269658427282, 24028.19815953374 12451.154294352575, 24058.96743142396 12389.615807871862, 24076.275018093882 12336.731229167028, 12247.308125834508, 24120.505746243296 12185.769639353795, 24128.198006985484 12107.885041575475, 24128.198006985484 12080.9620145369, 24123.390344021616 12058.846650462194, 24098.904966050724049.35187661439 11901.15416083817, 24024.3502920228 11865.57226023717, 23928.1980327943 11784.808144907995, 23892.621148464976 11743.4620145369, 23863.774941799937 11722.30813054967, 23806.082757351694 11663.654351573034, 23720.159539578257 11598.26956305985, 23666.659615872202 11587.692704539342, 23613.77486550599 11569.423356394811, 11608.8465936021276, 23502.236512539683 11615.57726417069, 23482.044213650523 11616.538911204381, 23464.736512539683 11581.923394541784, 23462.81347354139 11512.692590098424, 23460.8903821685923413.77482735900 11462.69266639237, 23322.42877328187 11426.154313426061, 23268.58271920472 11415.577340464635, 23237.813473541275 11415.577340464635, 23180.121148644976 11425.19266639237, 23145.6595615872202 11426.154313426061, 23132.813473541275 11415.577340464635, 23090.121148644976 11425.19266639237, 23059.121148644976 11425.19266639237, 23053.19793069154 11487.692628245397, 23024.242658840592 11474.23171946569, 22989.736550686655 11463.654313426061, 11448.26956305985, 22891.65953978257 11439.615884165807, 22858.00578439027 11439.615884165807, 22822.428849575816 11452.115807871862, 22810.890344021616 11464.615846018834, 22785.89038216859, 22766.659615872202 11567.90291209264, 22757.04417505355 11588.654351573034, 22746.467316983042 11635.76956305985, 22738.774941799937 11713.65427570908, 22745.50578439027 11742.50036750321, 21177.30813054967 22821.467316983042 11848.269524912877, 22853.198006985484 11884.000329356237, 22885.890344021616 11923.269753794713, 22899.35194761362 11948.269753794774, 22893.5826407638022842.621148644976 12067.500443797155, 22849.35187661439 12087.692742686315, 12113.654294352575, 22969.544251797495 12117.500481944127, 23023.390305874644 12130.000462807641, 23040.69782544566 12140.577321391149, 23090.69793069154 12181.9233908982702, 23116.65953957825, 23230.1595615872202 11567.90291209264, 23228.121148644976 11568.65953978257, 23206.89038216859, 23200.890344021616 11572.692590098424, 23188.774941799937 11582.30813054967, 23176.6595615872202 11583.846593241735, 23154.121148644976 11584.615884165807, 23132.242658840592 11585.23171946569, 23110.154313426061, 23114.73647439271 12098.269601206823, 23119.35176217347 12098.269601206823, 23119.928735134898 12092.5005200911, 23191.659654019175 12037.692704539342, 23205.12110318003 1, 23202.1595615872202 12037.115731577916, 23220.505746243296 12037.692704539342, 23207.04417550355 12037.692704539342, 23206.082642910776 12069.423508982702, 23218.58268105775 12096.346536021276, 12152.115826945348 23211.851838467417 12154.038834910436, 23193.58271920472 12149.23171946569), "geometryCompression": 2.25, "centroid": { "x": 23452.105097064356, "y": 11971.324686548254 }, "area": 997379.0, "areaUnit": "pixels", "perimeter": 5975.0, "perimeterUnit": "pixels", "term": [ "528044" ], "nbComments": 0, "cropURL": "https://research.cytomine.be/api/userannotation/143584607/crop.png", "smallCropURL": "https://research.cytomine.be/api/userannotation/143584607/crop.png?maxSize=256", "url": "https://research.cytomine.be/api/userannotation/143584607/crop.png", "imageURL": "https://research.cytomine.be/#/project/528050/image/528132/annotation/143584607", "reviewed": false }
```

Annotations in Cytomine REST API




<https://demo.cytomine.be/api/userannotation/143584607/crop.png>




(5556 x 4672 pixels)

<https://demo.cytomine.be/api/userannotation/143584607/crop.png?maxSize=512>




(512 x 430 pixels)

Annotations in Cytomine REST API



<https://demo.cytomine.be/api/userannotation/143584607/crop.png&draw=true>




<https://demo.cytomine.be/api/userannotation/143584607/crop.png?maxSize=512&mask=true>




Other API examples: get a tile (window)


<https://demo.cytomine.be/api/imageinstance/153042325/window-3000-3000-512-512.png?projection=max>



raw sub-image (window)



maximum projection of the window




(Support of hyperspectral images
is under development - Cytomine
ULiège R&D Edition)

Important note: In Cytomine coordinate system (0,0) is bottom left corner
 $(x,y) \rightarrow (x, \text{image.height}-y)$

Complete Cytomine HTTPS REST API Reference

- <https://doc.uliege.cytomine.org/dev-guide/api/reference>

In-house specification schema (move to OpenAPI ongoing)
 Still under validation (soon on doc.cytomine.org)



The screenshot shows a detailed API documentation page for the 'Get an annotation crop' endpoint. The page is structured as follows:

- Left Sidebar:** A navigation menu titled 'Data scientist Guide' containing links for 'API client for Python', 'API client for Java', 'API client for Javascript', and 'API reference'. The 'API reference' link is currently active.
- Header:** 'ANNOTATION' and 'Get an annotation crop'.
- Request Section:**
 - Path Parameters:** id (int64) - The annotation id; format (enum: jpg | png) - The expected image format.
 - Query-String Parameters:** type (string), draw (boolean), mask (boolean), alphaMask (boolean), maxSize (int32), zoom (int32), increaseArea (number), complete (boolean), colormap (string), inverse (int32), contrast (number), gamma (number), bits (int32).
- Response Section:** Shows status codes 200, 400, 401, 403, 404. Below this, there's an 'Image' example labeled 'binary' which is described as a 'Binary image'.
- Bottom Right:** A dropdown menu for 'image/jpeg' and a link for 'Single line description'.

Cytomine API clients (Python/Java/Javascript)

- **Goal:** to manipulate Cytomine resources as regular objects within your favorite language
 - Clients wrap the HTTP API endpoints:
 - Python: <https://github.com/cytomine/Cytomine-python-client>
 - Java: <https://github.com/cytomine/Cytomine-java-client>
 - JavaScript: <https://github.com/cytomine/Cytomine-js-client>
 - "Project" example with Python client:
 - List projects I have access to
 - `projects = ProjectCollection().fetch()`
 - Get a specific project (if I have rights)
 - `project = Project().fetch(80741)`
 - Update project info (if I have rights)
 - `project.name = "My new name"`
 - `project.update()`
 - Delete project (if I have rights)
 - `project.delete()`
 - Add new project (if I have rights)
 - `other_project = Project(name="My new project").save()`

Cytomine Python API client

- <https://doc.cytomine.org/dev-guide/clients/python/usage#models>

Data scientist Guide ▶

API client for Python ▾

Installation

Usage

Authentication

[Externalize credentials](#)

Models

Fetch a resource

Add a resource

Update a resource

Delete a resource

Fetch embedded resources

Collections

Filtering

Special case: Annotations

Search annotations

Multiple adds in one request

Verbose mode

Examples

Reference

API client for Java ▶

API client for Javascript ▶

API reference

Models

A Cytomine **resource** (such as a project *STUDY-01* or an image *cell01.tif*) is an instance of a **domain** (in the examples: a project, an image).

In the Cytomine API client for Python, a **resource is a Python object** which is an instance of a **model class** describing its domain. Each object (thus an instance of a model):

- has a set of attributes corresponding to resource attributes
- is managed itself through `fetch()`, `save()`, `update()` and `delete()` methods that communicate with the Cytomine server
- has some utilities to manage HTTP and JSON technical details

Examples of models are `Project`, `ImageInstance`, `Annotation`, `Term`, `User`, `CurrentUser`, etc.

Fetch a resource

To get data from server, we first need to fetch the resource. Indeed, when you create a new object (such as `myproject = Project()`), the object (here `myproject`) is empty. You need to populate this object with data coming from Cytomine.

TIP

The next sections show examples to fetch either a project or an image. The sections are very similar: the API client library is indeed very systematic.

Example: fetch a project

Supposing the Cytomine server has a project called *STUDY-01* with an identifier (ID) 42, the project resource is fetched in Python with:

```

1  from cytomine.models import Project
2
3  # ... Assume we are connected
4
5  myproject = Project().fetch(id=42)
6
7  # Will print STUDY-01
8  print(myproject.name)

```

py

Cytomine Java API client

- <https://doc.cytomine.org/dev-guide/clients/java/usage#models>

Data scientist Guide ▶

API client For Python ▶

API client for Java ▾

Installation

Usage

Authentication

Externalize credentials

Models

Fetch a resource

Add a resource

Update a resource

Delete a resource

Fetch embedded resources

Collections

Filtering

Special case: Annotations

Search annotations

API client for Javascript ▶

API reference

Models

A Cytomine **resource** (such as a project *STUDY-01* or an image *cell01.tif*) is an instance of a **domain** (in the examples: a project, an image).

In the Cytomine API client for Java, a **resource is a Java object** which is an instance of a **model class** describing its domain. Each object (thus an instance of a model):

- has a set of attributes corresponding to resource attributes
- is managed itself through `fetch()`, `save()`, `update()` and `delete()` methods that communicate with the Cytomine server
- has some utilities to manage HTTP and JSON technical details

Examples of models are `Project`, `ImageInstance`, `Annotation`, `Term`, `User`, `CurrentUser`, etc.

Fetch a resource

To get data from server, we first need to fetch the resource. Indeed, when you create a new object (such as `Project myproject = new Project()`), the object (here `myproject`) is empty. You need to populate this object with data coming from Cytomine.

TIP

The next sections show examples to fetch either a project or an image. The sections are very similar: the API client library is indeed very systematic.

Example: fetch a project

Supposing the Cytomine server has a project called *STUDY-01* with an identifier (ID) 42, the project resource is fetched in Java with:

Java	Groovy	java
<pre>1 import be.cytomine.client.models.Project; 2 3 // ... 4 // Assume we are connected 5 // ... 6 7 Long projectId = 42L; 8 Project myProject = new Project().fetch(projectId); 9 10 // Will print STUDY-01 11 String name = myProject.getStr("name"); 12 System.out.println(name);</pre>		

Behind the scene, the API client library has sent a `GET /api/project/42.json` request to the Cytomine server and populated `myproject` with the response content.

Only the project meta information have been retrieved (such as name, ontology identifier, settings,...). You need to perform other `fetch()` calls on other models to get the project images, annotations, project members, etc. See the [API reference](#) for all details.

Cytomine Java API client: example

- Java client examples:


<https://github.com/cytomine/Cytomine-java-client/tree/master/src/test/java/client>

```
Cytomine cytomite = new Cytomine(cytomiteCoreUrl, publickey, privatekey, "./");

println "List the images of your project:"
ImageInstanceCollection images = cytomite.getImageInstances(proj.getId());
for(int i=0;i<images.size();i++) {
    ImageInstance image = images.get(i)
    println image.getId()
}

println "List the annotations of your project:"
AnnotationCollection annotations = cytomite.getAnnotationsByProject(proj.getId());
for(int i=0;i<annotations.size();i++) {
    Annotation annotation = annotations.get(i)
    List termNameList = annotation.getList("term").collect{ termId ->
        return terms.list.find{it.id==termId}.get("name")
    }
    println annotation.getId() + " " + annotation.get("location") + " " + annotation.get("image") + " " + annotation.get("name")
}
```


Cytomine Java API client: e.g. access Cytomine data from another tool



<https://github.com/danyfel80/icytominne>

Cytomine Javascript API client: extend or create Web-UI components

- <https://doc.cytomine.org/dev-guide/clients/javascript/installation#how-to-install>
- Examples on <https://github.com/cytomine/Cytomine-js-client/tree/master/tests>



Single-Page Application, Written in VueJS, loaded in browser:

<https://github.com/cytomine/Cytomine-Web-UI>

Annotation	Users	Agreement
E01_Pronormoblast	masslaber, pkainz	67% (2)
E02_Basophilic Normoblast	masslaber, pkainz	100% (2)
E03_Polychromatophilic Normoblast	masslaber, pkainz	33% (1)
E04_Orthochromatic Normoblast	masslaber, pkainz	67% (2)
E05_Reticulocyte	masslaber, pkainz	33% (1)
Z01_Lymphocyte	masslaber, pkainz	67% (2)
Z02_Plasmacytoid Cell	masslaber, pkainz	67% (2)
Z03_Eosinophilic Granulocyte	masslaber, pkainz	67% (2)
Z04_Myeloblast	masslaber, pkainz	67% (2)
Z05_Neutrophilic Granulocyte	masslaber, pkainz	67% (2)
Z06_Normoblast	masslaber, pkainz	67% (2)
Z07_Proerythroblast	masslaber, pkainz	67% (2)
Z08_Reticular Cell	masslaber, pkainz	67% (2)
Z09_Unknown	masslaber, pkainz	67% (2)

Cytomine IRIS (Inter-observer Reliability Study module) by Ph. Kainz

<https://github.com/cytomine/Cytomine-IRIS> (not using JS client)

<https://www.youtube.com/watch?v=S0PhEJmqImA>

Cytomine Python client : quick demo

- Install Cytomine-Python-client or use “pre-built” Docker image
<https://doc.cytomine.org/dev-guide/clients/python/installation>
- Authentication (using host, public & private keys from [Account page](#))
- Fetch resources (e.g. project, ontology, terms, annotations, tile)
- Add resources (e.g. annotations, properties, upload image)

```
(base) maree@XPS-15-7590:~/cytome/neubias-webinar$ sudo docker run -it -v "$PWD":/usr/src/myapp -w /usr/src/myapp cytomeuliege/software-python3-base:latest python webinar-get-data.py
[sudo] password for maree:
[2021-05-11 10:38:49,089][INFO] [GET] [currentuser] CURRENT USER - 701 : rmaree | 200 OK
[currentuser] CURRENT USER - 701 : rmaree
701
None
-----
[2021-05-11 10:38:50,115][INFO] [GET] [project collection] 527 objects | 200 OK
[project collection] 527 objects
[2021-05-11 10:38:50,189][INFO] [GET] [project] 10688948 : CHALLENGE-CAMELYON16-TRAIN | 200 OK
Project NAME: CHALLENGE-CAMELYON16-TRAIN | Ontology: 10688937
-----
[2021-05-11 10:38:50,275][INFO] [ontology] 10688937 : CAMELYON16-TRAIN | 200 OK
[2021-05-11 10:38:50,359][INFO] [GET] [term collection] 5 objects | 200 OK
Term ID: 528906181 | Name: Nucleus
Term ID: 528905747 | Name: ROI
Term ID: 10694501 | Name: Tumor
Term ID: 119767342 | Name: Sample
Term ID: 120684669 | Name: Sample_pred
-----
[2021-05-11 10:38:51,116][INFO] [GET] [imageinstance collection] 271 objects | 200 OK
Name: tumor_111.tif | Width: 97792 | Height: 220672 | ImageID: 125220817
Name: Tumor_009.tif | Width: 97792 | Height: 217088 | ImageID: 27712918
Name: Tumor_008.tif | Width: 97792 | Height: 221184 | ImageID: 27712752
Name: Tumor_007.tif | Width: 97792 | Height: 221184 | ImageID: 27712586
Name: Tumor_006.tif | Width: 97792 | Height: 221184 | ImageID: 27712408
Name: Tumor_005.tif | Width: 97792 | Height: 219648 | ImageID: 27712229
Name: tumor_110.tif | Width: 94208 | Height: 71680 | ImageID: 27712042
Name: tumor_109.tif | Width: 131072 | Height: 89600 | ImageID: 27711876
Name: tumor_108.tif | Width: 126976 | Height: 107520 | ImageID: 27711710
```

Cytomine Python client : applying a pre-trained Stardist model

- Stardist + Cytomine Python client to communicate with Cytomine server

1. Get ROIs (images with alphamask) from Cytomine server, save local PNG file

```

from csbdeep.utils import Path, normalize
from stardist import random_label_cmap
from stardist.models import StarDist2D
from cytomine import cytomine, models, CytomineJob
from cytomine.models import Annotation, AnnotationTerm, AnnotationCollection, ImageInstanceCollection, Job
from PIL import Image

roi_annotations = AnnotationCollection()
roi_annotations.project = conn.parameters.cytomine_id_project
roi_annotations.term = conn.parameters.cytomine_id_roi_term
roi_annotations.image = id_image
roi_annotations.showWKT = True
roi_annotations.fetch()
print(roi_annotations)
#Go over ROI in this image
for roi in roi_annotations:
    #Get Cytomine ROI coordinates for remapping to whole-slide
    #Cytomine cartesian coordinate system, (0,0) is bottom left corner
    roi_geometry = wkt.loads(roi.location)
    minx=roi_geometry.bounds[0]
    miny=roi_geometry.bounds[3]
    #Dump ROI image into local PNG file
    roi_path=os.path.join(working_path,str(roi_annotations.project)+'/'+str(roi_annotations.image)+'/'+str(roi.id))
    roi_png_filename=os.path.join(roi_path+'/'+str(roi.id)+'.png')
    roi.dump(dest_pattern=roi_png_filename,mask=True,alpha=True)

```

Cytomine Python client : applying a pre-trained Stardist model

- Stardist + Cytomine Python client to communicate with Cytomine server

2. Apply pre-trained Stardist model, convert & upload annotations to Cytomine

```

labels, details = model.predict_instances(img,
                                         prob_thresh=conn.parameters.stardist_prob_t,
                                         nms_thresh=conn.parameters.stardist_nms_t)
print("Number of detected polygons: %d" %len(details['coord']))
cytominе_annotations = AnnotationCollection()
#Go over detections in this ROI, convert and upload to Cytomine
for pos,polygroup in enumerate(details['coord'],start=1):
    #Converting to Shapely annotation
    points = list()
    for i in range(len(polygroup[0])):
        #Cytomine cartesian coordinate system, (0,0) is bottom left corner
        #Mapping Stardist polygon detection coordinates to Cytomine ROI in whole slide image
        p = Point(minx+polygroup[1][i],miny-polygroup[0][i])
        points.append(p)

    annotation = Polygon(points)
    #Append to Annotation collection
    cytominе_annotations.append(Annotation(location=annotation.wkt,
                                           id_image=id_image,
                                           id_project=conn.parameters.cytominе_id_project,
                                           id_terms=[conn.parameters.cytominе_id_cell_term]))
#Send Annotation Collection (for this ROI) to Cytomine server in one http request
ca = cytominе_annotations.save()

```

Cytomine Python client : applying a pre-trained Stardist model

- Stardist + Cytomine Python client to communicate with Cytomine server
 - 2. Apply pre-trained Stardist model, convert & upload annotations to Cytomine


OK but...

- It runs with user credentials and create data in [Cytomine user own layer \(UserAnnotation\)](#)
 - It is neither traceable nor reproducible
 - Executed on my own computer with specific OS/libraries
 - Code is not versioned
 - Parameter values are not saved
- We need more formalism to improve usability & reproducibility

```
#Send Annotation Collection (for this ROI) to Cytomine server in one http request
ca = cytomine_annotations.save()
    id_project=conn.parameters.cytomine_id_project,
    id_terms=[conn.parameters.cytomine_id_cell_term]))
```

Cytomine “Apps” concepts: Software, Job, Parameters

<https://doc.cytomine.org/dev-guide/algorithms/>



Cytomine “Apps”: Parameter Descriptor and Containerized Execution Environment

e.g. https://github.com/cytomine/S_CellDetect_Stardist_HE_ROI

The diagram illustrates the components of a Cytomine “App”:

- Dockerfile:** A file defining the container environment. It specifies the base image as `cytominelab/software-python3-base:v2.2.0` and installs dependencies like TensorFlow and Stardist. It also adds configuration files and models to the `/models` directory.
- Descriptor:** A JSON file defining parameters. It includes entries for `stardist_prob_t` (Stardist Probability Threshold) and `stardist_nms_t` (Stardist Non-Maximum Suppression Overlap threshold).
- Code (e.g. Python):** A snippet of Python code for performing cell detection. It uses the `predict_instances` method from a model, processes results, and converts them into Cytomine ROI annotations.
- GitHub Repository:** A screenshot of the GitHub repository page for `S_CellDetect_Stardist_HE_ROI`. It shows the repository structure, including `Dockerfile`, `LICENSE`, `README.md`, `descriptor.json`, `run.py`, and several image files. A pull request titled “geektoroise integration to cytomite repository” is shown, along with a commit history and a screenshot of the repository interface.
- README.md:** A file containing the title `S_CellDetect_Stardist_ROI`.

Cytomine “Apps” Parameter Descriptor in JSON

e.g. [https://github.com/cytomine/S_Celldetect_Stardist_HE_ROI/blob/master\(descriptor.json](https://github.com/cytomine/S_Celldetect_Stardist_HE_ROI/blob/master(descriptor.json)

follows: <https://doc.cytomine.org/dev-guide/algorithms/descriptor-reference>

Execution command line + mandatory Cytomine parameters (host, private/public key, project id, software id)

+ algorithm specific parameters (e.g. thresholds,...)

```
{
  "name": "Celldetect_Stardist_HE_ROI",
  "container-image": {
    "image": "cytomine/s_celldetect_stardist_he_roi",
    "type": "singularity"
  },
  "schema-version": "cytominne-0.1",
  "description": "Cell (nuclei) detection using StarDist with versatile H&E pre-trained model",
  "command-line": "python run.py CYTOMINE_HOST CYTOMINE_PUBLIC_KEY CYTOMINE_PRIVATE_KEY CYTOMINE_ID_PROJECT CYTOMINE_ID_SOFTWARE CYTOMINE_ID_IMAGES
  "inputs": [
    {
      "id": "cytominne_host",
      "value-key": "@ID",
      "command-line-flag": "--@id",
      "name": "Cytomine host",
      "set-by-server": true,
      "optional": false,
      "type": "String"
    },
    {
      "id": "cytominne_public_key",
      "value-key": "@ID",
      "command-line-flag": "--@id",
      "name": "Cytomine public key",
      "set-by-server": true,
      "optional": false,
      "type": "String"
    },
    {
      "id": "cytominne_private_key",
      "value-key": "@ID",
      "command-line-flag": "--@id",
      "name": "Cytomine private key",
      "set-by-server": true,
      "optional": false,
      "type": "String"
    },
    {
      "id": "cytominne_id_project",
      "value-key": "@ID",
      "command-line-flag": "--@id",
      "name": "Cytomine project ID",
      "set-by-server": true,
      "optional": false
    }
  ],
  "outputs": [
    {
      "id": "stardist_prob_t",
      "value-key": "@ID",
      "command-line-flag": "--@id",
      "name": "StarDist Probability Threshold",
      "description": "Probability Threshold in range [0.0, 1.0] - higher values lead to fewer segmented objects, but will likely result in more false positives",
      "default-value": 0.5,
      "set-by-server": false,
      "optional": true,
      "type": "Number"
    },
    {
      "id": "stardist_nms_t",
      "value-key": "@ID",
      "command-line-flag": "--@id",
      "name": "StarDist Non-Maximum Suppression Overlap threshold",
      "description": "Overlap Threshold in range [0.0, 1.0] - higher values allow segmented objects to overlap substantially.",
      "default-value": 0.5,
      "set-by-server": false,
      "optional": true,
      "type": "Number"
    }
  ]
}
```

The JSON descriptor is parsed by Cytomine server to create a Software in Cytomine database tables
 → “App” resource: <https://demo.cytomine.be/api/software/155177168.json>
 → “App” webpage: <https://demo.cytomine.be/#/algorithm/155177168>

See <https://doc.cytomine.org/dev-guide/algorithms/descriptor-reference#complete-example>
<https://doc.cytomine.org/dev-guide/algorithms/write-app>

Cytomine “Apps” Environment Descriptor (Dockerfile)

e.g. https://github.com/cytomine/S_CellDetect_Stardist_HE_ROI/blob/master/Dockerfile

```
FROM cytomine/software-python3-base:v2.2.0
#
# -----
# Install Stardist and tensorflow
RUN pip install tensorflow==2.2.0
RUN pip install stardist==0.6.0
RUN mkdir -p /models && \
    cd /models && \
    mkdir -p 2D_versatile_HE
ADD config.json /models/2D_versatile_HE/config.json
ADD thresholds.json /models/2D_versatile_HE/thresholds.json
ADD weights_best.h5 /models/2D_versatile_HE/weights_best.h5
RUN chmod 444 /models/2D_versatile_HE/config.json
RUN chmod 444 /models/2D_versatile_HE/thresholds.json
RUN chmod 444 /models/2D_versatile_HE/weights_best.h5

#
# -----
# Install scripts
ADD descriptor.json /app/descriptor.json
RUN mkdir -p /app
ADD run.py /app/run.py

ENTRYPOINT ["python3", "/app/run.py"]
```

Base container image that includes Cytomine Python client

Specific libraries & files (e.g. pre-trained model) needed by your code and installed in the container image

Command to be executed when container image is instanciated

The Dockerfile is used by Docker to build a reproducible image of your software environment.
It allows to execute almost anything e.g. codes from:



Cytomine “Apps” versioning

- Versioning is achieved by creating Github release tags for important new code versions
- Automatic Builds configuration (DockerHub ↔ Github) are used to create Docker images for each new release and automatically import them into Cytomine server (“trusted sources” configuration)

The image shows a GitHub repository interface and a Cytomine application management interface side-by-side.

GitHub Repository:

- Repository: Cytomine-ULiege / S_CellDetect_Stardist_HE_ROI
- Tags tab selected.
- Tags listed:
 - v1.0.3 (released on May 12, 2020, commit 01e241d, zip, tar.gz, Notes)
 - v1.0.2 (released on May 11, 2020, commit 10f1116, zip, tar.gz, Notes)
 - v1.0.1 (released on May 11, 2020, commit 7cefadb, zip, tar.gz, Notes)

Cytomine Application List:

- General tab selected.
- Algorithms tab selected.
- Search bar: Stats-TernArea (v1.1.1)
- Table of applications:

Name	Version	Runnable	Status
Stats-TernArea (v1.1.1)	Last release	Yes	Enabled
CellDetect_Stardist_HE_ROI (v1.0.3)	Last release	Yes	Enabled
CellDetect_Stardist_HE_ROI (v1.0.2)	Deprecated	Yes	Enabled
CellDetect_Stardist_HE_ROI	No information	No	Enabled
Segment-ML-UNet-Binary-Pred	No information	No	Disabled
Segment-ML-UNet-Binary-Train (v0.1.1)	Last release	Yes	Disabled
Segment-ML-UNet-Binary-Train (v0.1.0)	Deprecated	Yes	Disabled
Segment-ML-UNet-Binary-Train	No information	No	Disabled
S_Zebrafish_Head_Operculum_UNet_Segmentation (v3.0.4)	Last release	Yes	Disabled
S_Zebrafish_Head_Operculum_UNet_Segmentation (v3.0.3)	Deprecated	Yes	Disabled
S_Zebrafish_Head_Operculum_UNet_Segmentation (v3.0.2)	Deprecated	Yes	Disabled

Cytomine “Apps” execution logging

- Traceability is achieved by logging software stack trace

Execution log

Hide

```
[2020-12-03 16:25:40,460][INFO] [GET] [currentuser] CURRENT USER - 525457679 : JOB[rmree], 2020-12-03 16:25:30:631 | 200 OK
[2020-12-03 16:25:40,525][INFO] [GET] [project] 138124787 : DEMONSTRATION-VARIOUS | 200 OK
[2020-12-03 16:25:40,594][INFO] [GET] [software] 155177168 : CellDetect Stardist_HE_ROI | 200 OK
[2020-12-03 16:25:40,642][INFO] [GET] [softwareparameter collection] 12 objects | 200 OK
[2020-12-03 16:25:40,712][INFO] [GET] [job] 525457631 : None | 200 OK
[2020-12-03 16:25:40,712][INFO] Job (id:525457631) status update: "None" (status: RUNNING, progress: 0%)
[2020-12-03 16:25:40,839][INFO] [PUT] [job] 525457631 : None | 200 OK
[2020-12-03 16:25:40,840][INFO] Job (id:525457631) status update: "Initialization..." (status: RUNNING, progress: 0%)
[2020-12-03 16:25:40,960][INFO] [PUT] [job] 525457631 : None | 200 OK
Using TensorFlow backend.
WARNING:tensorflow:From /usr/local/lib/python3.7/site-packages/tensorflow_core/python/ops/resource_variable_ops.py:1630: calling BaseResourceVariable.__init__ (from tensorflow.python.ops.resource_variable_ops) with constraint is deprecated and will be removed in a future version.
Instructions for updating:
If using Keras pass *_constraint arguments to layers.
[2020-12-03 16:25:41,795][WARNING] From /usr/local/lib/python3.7/site-packages/tensorflow_core/python/ops/resource_variable_ops.py:1630: calling BaseResourceVariable.__init__ (from tensorflow.python.ops.resource_variable_ops) with constraint is deprecated and will be removed in a future version.
Instructions for updating:
If using Keras pass *_constraint arguments to layers.
WARNING:tensorflow:From /usr/local/lib/python3.7/site-packages/keras/backend/tensorflow_backend.py:4070: The name tf.nn.max_pool is deprecated. Please use tf.nn.max_pool2d instead.

[2020-12-03 16:25:41,839][WARNING] From /usr/local/lib/python3.7/site-packages/keras/backend/tensorflow_backend.py:4070: The name tf.nn.max_pool is deprecated. Please use tf.nn.max_pool2d instead.

2020-12-03 16:25:42.396999: W tensorflow/stream_executor/platform/default/dso_loader.cc:55] Could not load dynamic library 'libcuda.so.1'; dlerror: libcuda.so.1: cannot open shared object file: No such file or directory; LD_LIBRARY_PATH: ./singularity.d/libs
2020-12-03 16:25:42.397090: E tensorflow/stream_executor/cuda/cuda_driver.cc:318] failed call to cuInit: UNKNOWN ERROR (303)
2020-12-03 16:25:42.397141: I tensorflow/stream_executor/cuda/cuda_diagnostics.cc:156] kernel driver does not appear to be running on this host (cytominer-slurm): /proc/driver/nvidia/version does not exist
2020-12-03 16:25:42.397453: I tensorflow/core/platform/cpu_feature_guard.cc:142] Your CPU supports instructions that this TensorFlow binary was not compiled to use: FMA
2020-12-03 16:25:42.458346: I tensorflow/core/platform/profile_utils/cpu_utils.cc:94] CPU Frequency: 2600085000 Hz
2020-12-03 16:25:42.469674: I tensorflow/compiler/xla/service/service.cc:168] XLA service 0x55c4cbf43880 initialized for platform Host (this does not guarantee that XLA will be used). Devices:
2020-12-03 16:25:42.469876: I tensorflow/compiler/xla/service/service.cc:176] StreamExecutor device (0): Host, Default Version
Loading network weights from 'weights_best.h5'.
Loading thresholds from 'thresholds.json'.
Using default values: prob_thresh=0.692478, nms_thresh=0.3.
[2020-12-03 16:25:42,993][INFO] [imageinstance collection] 5 objects | 200 OK
[2020-12-03 16:25:42,994][INFO] Job (id:525457631) status update: "Running detection on image (1/1)." (status: RUNNING, progress: 0%)
[2020-12-03 16:25:43,108][INFO] [PUT] [job] 525457631 : None | 200 OK
[2020-12-03 16:25:43,189][INFO] [GET] [annotation collection] 1 objects | 200 OK
[annotation collection] 1 objects
-----ROI-----
ROI Geometry from Shapely: POLYGON ((3.341142933359549 508.4867133364421, 675.7250540661721 508.4867133364421, 675.7250540661721 4.675006793473358, 3.341142933359549 4.675006793473358, 3.341142933359549 508.4867133364421))
ROI Bounds
(3.341142933359549, 4.675006793473358, 675.7250540661721, 508.4867133364421)
roi_png_filename: /home/cytomine/525457631/138124787/525455833/525455855/525455855.png
[2020-12-03 16:25:44,041][INFO] File downloaded successfully from https://research.cytomine.be/api/userannotation/525455855/alphamask.png
WARNING:tensorflow:From /usr/local/lib/python3.7/site-packages/keras/backend/tensorflow_backend.py:422: The name tf.global_variables is deprecated. Please use tf.compat.v1.global_variables instead.

Normalizing image channels independently.
----- Processing ROI file 0: /home/cytomine/525457631/138124787/525455833/525455855/525455855.tif
[2020-12-03 16:25:44,215][WARNING] From /usr/local/lib/python3.7/site-packages/keras/backend/tensorflow_backend.py:422: The name tf.global_variables is deprecated. Please use tf.compat.v1.global_variables instead.

Number of detected polygons: 60
-----[2020-12-03 16:25:45,211][INFO] [POST] algoannotations added | 200 OK
[2020-12-03 16:25:47,052][INFO] [POST] algoannotations 525457740,52545766,525457784,52545788,525457932,525457995,525458061,525458114,525458141,525458170,525458213,525458252,525458291,525458357 added | 200 OK
[2020-12-03 16:25:47,217][INFO] [POST] algoannotations 525457756,525457810,525457907,525457955,525458023,525458088,525458133,525458187,525458226,525458285,525458319,525458374,525458412,525458450 added | 200 OK
[2020-12-03 16:25:47,281][INFO] [POST] algoannotations 525457757,525457794,525457853,525457919,525457971,525458006,525458049,525458105,525458162,525458244,525458309,525458349,525458387,525458438,525458473 added | 200 OK
[2020-12-03 16:25:47,282][INFO] [POST] algoannotations 525457758,525457824,525457875,525457981,525458037,525458075,525458135,525458208,525458269,525458334,525458401,525458462,525458495,525458508 added | 200 OK
[2020-12-03 16:25:47,429][INFO] Job (id:525457631) status update: "Finished." (status: SUCCESS/TERMINATED, progress: 100%)
[2020-12-03 16:25:47,556][INFO] [PUT] [job] 525457631 : None | 200 OK
[2020-12-03 16:25:47,557][INFO] Job (id:525457631) status update: "Job successfully terminated" (status: SUCCESS/TERMINATED, progress: 100%)
[2020-12-03 16:25:47,675][INFO] [PUT] [job] 525457631 : None | 200 OK
```

Cytomine “Apps” reproducibility

- Reproducibility is further achieved by storing each execution parameter values (done by Cytomine-Core in its internal database for each execution)

Segment-ML-MultiOutputET-Pred-BI (v0.2.7)		4	rhoyoux	Mar 29, 2021 10:51 AM	Success
Status comment	Job successfully terminated				
Execution duration	18 minutes				
Parameters	Hide	Name	Value	Type	
Zoom level		5		Number	
Tile size		256		Number	
Tile overlap		0		Number	
Number of jobs		8		Number	
Union enabled		true		Boolean	
Min annotations area		0		Number	
Max annotations area		-1		Number	
Minimum standard deviation for tile filtering		0		Number	
Maximum mean for tile filtering		255		Number	
Training job identifier		207885		Cytomine domain	
Prediction step		1		Number	
Cytomine host				String	
Cytomine public key				String	
Cytomine private key				String	
Cytomine software id		202996		Number	
Cytomine project id		1397		Number	

→ Reproducible execution of Cytomine “Apps”

- From Cytomine WebUI

The screenshot shows the 'Launch new analysis' dialog box. At the top left is a 'Analysis' button with a three-bars icon. To its right is a blue 'Launch new analysis' button. The main area contains the following information:

Launch new analysis

Algorithm: CellDetect_Stardist_HE_ROI (v1.0.4)

Name	Value
Cytomine Image IDs	CMU-1.tiff
Cytomine ROI term ID	ROI
Cytomine Cell term ID	Nucleus

Pre-filled parameters Hide

Stardist Probability Threshold	0.5
Stardist Non-Maximum Suppression Overlap threshold	0.5
Stardist Image Normalization Percentile Low	1
Stardist Image Normalization Percentile High	99.8

Cancel Launch new analysis



When you execute a Cytomine “App” from Cytomine Web-UI, Cytomine server sends parameter values to the “App” software environment (Docker container image) that is running on a “Processing Server”.
This “App” interacts with Cytomine using its embedded Cytomine Python client and upload results to the server

→ Reproducible execution of Cytomine “Apps”

From remote command-line (on your own computer/server/cluster):

- sudo docker pull cytomite/s_celldetect_stardist_he_roi:v1.0.4 ← Github/Dockerhub release tag
- sudo docker run -it cytomite/s_celldetect_stardist_he_roi:v1.0.4
--cytomin`e_host` <https://demo.cytomine.be> \
--cytomin`e_public_key` XXXXX \
--cytomin`e_private_key` XXXXX \
--cytomin`e_id_project` 6870 \
--cytomin`e_id_software` 2762 \
--cytomin`e_id_images` 6967 \
--cytomin`e_id_roi_term` 6880 \
--cytomin`e_id_cell_term` 6888 \
--stardist`_prob_t` 0.8 \
← Cytomine credentials from Account page
← Cytomine resource identifiers
← Your algorithm specific parameters for this execution

→ Then results are also accessible in Cytomine WebUI:
<https://demo.cytomine.be/#/project/10688948/analysis>



Cytomine and BIAFLOWS for data and computer scientists

Cytomine might look “complex” compared to regular desktop tools. But it allows to make image-based research truly **collaborative and reproducible** by sharing almost anything (images, annotations, apps, code, results, ...) **over the web**. Both end-users (WebUI) and computer scientists (API/clients) are able to manipulate, share, and publish all data.

Step-by-step documentation: <https://doc.cytomine.org/>

Any user issue ? <https://forum.image.sc/tag/cytomine>

Any code issue ? <https://github.com/cytomine/>



How to contribute ?

- Integrate your own “apps” and share them (thousands of potential users)
- Other ways: <https://doc.cytomine.org/community/how-to-contribute>

BIAFLOWS is a specific instance of Cytomine:

It uses same architecture/technologies, it comes with a slightly different UI, pre-populated workflows & datasets (2D+c+z+t); standardized inputs/outputs & metrics for benchmarking.



Acknowledgments

ULiège & Cytomine.coop : Ulysse Rubens, Romain Mormont, Navdeep Kumar, Remy Vandaele, Renaud Hoyoux, Grégoire Vincke, Pierre Geurts [cytomine.org]

Neubias and Comulis EU network : Sebastien Tosi, Benjamin Pavie, Volker Backer, Natasa Sladoje [neubias.org & comulis.eu]

Fundings (since 2010):

